

Troubleshooting (Taken from Bozzolla & Russel, 1999, Wagner)

Knife cuts every other section

- a. The advance has been set below the capabilities of the cutting edge. Increase the advance until serial sections are cut or use a sharper knife.

Failure to cut any sections

- a. Cantilever arm at end of fine advance.
- b. Dull knife.
- c. Block too soft.
- d. Knife or block not secure.
- e. Negative clearance angle.
- f. Wet block face.
- g. Vibrations.
- h. Temperature fluctuations.

Thickness variations from one entire section to the next

- a. Dull knife.
- b. Bumping of microtome.
- c. Drafts or temperature variations.
- d. Knife or block not secure.
- e. Block face too large or too soft.
- f. Wrong cutting speed.

Wrinkled sections

- a. Block face too large or too soft.
- b. Dirty or dull knife.
- c. Clearance angle too great.
- d. Water level too low.
- e. Cutting speed too fast.
- f. Knife not secure.

Compressed sections

- a. Block too soft.
- b. Cutting speed too fast.
- c. Inadequate expansion, try using.

Chatter

- a. High-frequency vibrations during sectioning, try a different cutting speed or clearance angle.
- b. Block too tall with small base.
- c. Dull knife or soft block.
- d. Block or knife not secure.

Specimen block lifts sections on return stroke

- a. Water level too high.
- b. Block face dirty, wet or hydrophilic. If block face is wet, wick dry with wedge of filter paper without touching knife edge.
- c. Clearance angle too small.
- d. Dirty knife or back of knife is wet.
- e. Static electricity on block face.

Block gets wet

- a. See a – e above.
- b. Block face too large.
- c. Cutting speed too low.

Sections fragged over knife edge

- a. Cutting speed too low.
- b. Water level too high.
- c. Clearance angle too low.
- d. Block too soft or a ragged edge of trapezoid preventing detachment.

Sections have holes

- a. Bubbles in resin.
- b. Incomplete infiltration with resin.
- c. Hard objects in specimen.

Specimen falls out of block

- a. Poor infiltration.
- b. Block too soft.

Sections have striations perpendicular to the knife

- a. Nick in knife edge. Move to a different region of knife edge or change knife.
- b. Dirt on knife edge.
- c. Knife damaged by hard region in specimen. Trim block to avoid hard region.

Sections do not form ribbons

- a. Top and bottom of trapezoid not parallel. Re-trim block.
- b. Water level wrong.
- c. Cutting speed too slow.
- d. Static electricity on block face.

Ribbon of sections curved

- a. Top and bottom of trapezoid not parallel. Try re-trimming block.
- b. Compression on one side of section.

Sections stick to eyelash probe

- a. Dirty eyelash probe.
- b. Bearing down on sections too much with eyelash probe.

Knife does not wet

- a. Add a drop of dilute Tween 20 solution to the boat.
- b. Use saliva to wet knife edge.

Sections hard to see

- a. Water level wrong.
- b. Illumination wrong.

Sections hard to move in boat

- a. Contamination in boat water. Change water.

Sections move away from grid

- a. Dirty grid.

Perpendicular regions with varied interference colors in sections

- a. Cutting edge not equally sharp across knife edge. Use different part of knife edge or change knife.

Irregular variations in interference colors throughout sections

- a. Uneven consistency between specimen and embedding material or within different regions of the specimen. Try to re-trim to include only areas with an even consistency.

Color variations occur in bands parallel to knife edge

- a. Low frequency vibrations.
- b. Knife or specimen not secure.
- c. Cutting speed too fast.
- d. Trapezoid needs to be re-trimmed.